

Detergents & Cleaners



Alkaline Detergents for Automatic Washers

Name	Part #	Description
Mild Alkaline Detergent	FHC-110CK-1	Instrument detergent for automatic washers. 1 gallon plastic bottle packed 4 per case.
Mild Alkaline Detergent	FHC-110CK-5	Instrument detergent for automatic washers. 5 gallon plastic jug w/ handle.
Mild Alkaline Detergent	FHC-110CK-15	Instrument detergent for automatic washers. 15 gallon plastic drum.
Mild Alkaline Detergent	FHC-110CK-30	Instrument detergent for automatic washers. 30 gallon plastic drum.
Mild Alkaline Detergent	FHC-110CK-55	Instrument detergent for automatic washers. 55 gallon plastic drum.
Mild Alkaline Detergent - Super Concentrated	FHC-110SCLF-1	Super concentrated detergent for automatic instrument washers. 1 gallon plastic bottle packed 4 per case.
Mild Alkaline Detergent - Super Concentrated	FHC-110SCLF-2.5	Super concentrated detergent for automatic instrument washers. 2.5 gallon plastic bottles packed by case.
Mild Alkaline Detergent - Super Concentrated	FHC-110SCLF-5	Super concentrated detergent for automatic instrument washers. 5 gallon plastic jug with handle.
Mild Alkaline Detergent - Super Concentrated	FHC-110SCLF-55	Super concentrated detergent for automatic instrument washers. 55 gallon plastic drum.

Detergents & Cleaners



Neutralizer for Automatic Washers

Name	Part #	Description
Acidic Neutralizer	FHC-200-5	Neutralizes residual alkaline detergent in the rinse cycle. 5 gallon plastic jug.
Acidic Neutralizer	FHC-200-15	Neutralizes residual alkaline detergent in the rinse cycle. 15 gallon plastic drum.
Acidic Neutralizer	FHC-200-30	Neutralizes residual alkaline detergent in the rinse cycle. 30 gallon plastic drum.
Acidic Neutralizer	FHC-200-55	Neutralizes residual alkaline detergent in the rinse cycle. 55 gallon plastic drum.
Acidic Neutralizer - Super Concentrated	FHC-200SCLF-2.5	Neutralizes residual alkaline detergent in the rinse cycle. 2 X 2.5 gallon plastic jugs.
Acidic Neutralizer - Super Concentrated	FHC-200SCLF-5	Neutralizes residual alkaline detergent in the rinse cycle. 5 gallon plastic jug.
Acidic Neutralizer - Super Concentrated	FHC-200SCLF-15	Neutralizes residual alkaline detergent in the rinse cycle. 15 gallon plastic drum.
Acidic Neutralizer - Super Concentrated	FHC-200SCLF-30	Neutralizes residual alkaline detergent in the rinse cycle. 30 gallon plastic drum.

Detergents & Cleaners



Lubricant for Automatic Washers

Name	Part #	Description
Lubricant Concentrate	FHC-300-5	Lubricant to protect surgical instruments from corrosion or staining in lubrication cycle. 5 gallon plastic jug.
Lubricant Concentrate	FHC-300-15	Lubricant to protect surgical instruments from corrosion or staining in lubrication cycle. 15 gallon plastic jug.
Lubricant Super Concentrate	FHC-300SCLF-2.5	Lubricant to protect surgical instruments from corrosion or staining in lubrication cycle. 2 x 2.5 gallon plastic jugs per box.
Lubricant Super Concentrate	FHC-300SCLF-5	Lubricant to protect surgical instruments from corrosion or staining in lubrication cycle. 5 gallon plastic jug.
Lubricant Super Concentrate	FHC-300SCLF-15	Lubricant to protect surgical instruments from corrosion/ staining in lubrication cycle. 15 gallon plastic jug.
Instrument Lubricant Milk	FHC-301-1	Lubricant for manual and machine use. Lubricates and protects instruments. 1 gallon plastic bottle/4 per case.
Instrument Lubricant Milk	FHC-301-5	Lubricant for manual and machine use. Lubricates and protects instruments. 5 gallon HDPE plastic jug w/ handle.
Instrument Lubricant Milk	FHC-301-15	Lubricant for manual and machine use. Lubricates and protects instruments. 15 gallon HDPE plastic drum.
Instrument Lubricant Milk	FHC-301-30	Lubricant for manual and machine use. Lubricates and protects instruments. 30 gallon HDPE plastic drum.
Instrument Lubricant Milk	FHC-301-55	Lubricant for manual and machine use. Lubricates and protects instruments. 55 gallon HDPE plastic drum.
Instrument Spray Lubricant	FHC-302S-8	A convenient, ready to use spray lubricant to protect and extend instrument life. 8 fl. oz. HDPE spray bottles / 12 per case.
Instrument Spray Lubricant	FHC-302S-16	A convenient, ready to use spray lubricant to protect and extend instrument life. 16 fl. oz. HDPE spray bottles / 8 per case.

Detergents & Cleaners



Neutral pH Detergent for Manual & Automatic Cleaning

Name	Part #	Description
Manual Neutral pH Detergent	FHC-401-1	Neutral pH detergent for pre-soaking and manual cleaning of instruments and glassware. 1 gallon plastic bottle/ 4 per case including pump dispenser.
Manual Neutral pH Detergent	FHC-401-5	Neutral pH detergent for pre-soaking and manual cleaning of instruments and glassware. 5 gallon plastic jug w/ handle.
Neutral pH Detergent	FHC400-1	Neutral pH low foam detergent for automatic and manual cleaning. 1 gallon plastic bottle/ 4 per case.
Neutral pH Detergent	FHC400-5	Neutral pH low foam detergent for automatic and manual cleaning. 5 gallon plastic jug w/ handle.
Neutral pH Detergent	FHC400-15	Neutral pH low foam detergent for automatic and manual cleaning. 15 gallon plastic drum.
Neutral pH Detergent	FHC400-30	Neutral pH low foam detergent for automatic and manual cleaning. 30 gallon plastic drum.
Neutral pH Power Detergent	FHC400-P5	Neutral pH powder detergent for automatic and manual cleaning. 5 lbs. pails/ 6 per case.
Neutral pH Power Detergent	FHC400-P25	Neutral pH powder detergent for automatic and manual cleaning. 25 lbs. fiber drum
Neutral pH Power Detergent	FHC400-P85	Neutral pH powder detergent for automatic and manual cleaning. 30 gallon plastic drum.
Neutral pH Cleaner - Super Concentrated	FHC400-2.5	Neutral pH super concentrated low foam cleaner for manual and automatic cleaning. 2 x 2.5 gallon plastic bottles per case.
Neutral pH Cleaner - Super Concentrated	FHC400-5	Neutral pH super concentrated low foam cleaner for manual and automatic cleaning. 5 gallon plastic jug w/ handle.
Neutral pH Cleaner - Super Concentrated	FHC400-15	Neutral pH super concentrated low foam cleaner for manual and automatic cleaning. 15 gallon plastic jug w/ handle.
Neutral pH Cleaner - Super Concentrated	FHC400-30	Neutral pH super concentrated low foam cleaner for manual and automatic cleaning. 30 gallon plastic drum.

Detergents & Cleaners



Neutral pH Pre-Cleaner & Enzymatic Detergent

Name	Part #	Description
Neutral pH Pre-Cleaner	FHC-402S-22	Enzymatic foam-spray neutral pH pre-cleaner guards against post-op residual bio-burden build up. 22 fl. oz. sprayer bottle & sprayer/ 12 per case.
Neutral pH Pre-Cleaner	FHC-402S-Gal	Enzymatic foam-spray neutral pH pre-cleaner guards against post-op residual bio-burden build up.. 4 x 1 gallon bottle packed 4 per case. (Refill for 22 oz.)
Enzymatic Pre-Cleaner	FHC-402S-5	Enzymatic foam-squeeze-spraypre-cleaner guards against post-op residual bio-burden build up. 48 x 5oz bottles in a case.
Neutral pH Enzymatic Detergent	FHC-404BB-1	Surg-ENZ® Bag-in a Box reduces cleaning time by rapid, efficient removal of stubborn organic encrustation on instruments and scope washers. 2 x 1 gal per case.
Neutral pH Enzymatic Detergent - Regular Foam	FHC-404C-1	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 4 x 1 gallon bottles per case
Neutral pH Enzymatic Detergent - Regular Foam	FHC-404C-5	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 5 gallon plastic jug w/ handle.
Neutral pH Enzymatic Detergent - Regular Foam	FHC-404C-15	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 15 Gallon plastic drum.
Neutral pH Enzymatic Detergent - Regular Foam	FHC-404C-30	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 30 gallon plastic drum.
Neutral pH Enzymatic Detergent - Regular Foam	FHC-404C-55	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 55 gallon plastic drum.

Detergents & Cleaners



Neutral pH Enzymatic Detergent (continued)

Name	Part #	Description
Neutral pH Enzymatic Detergent - Low Foam	FHC-404LF-1	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation on instruments and scope washers. 4 x 1 gallon bottles per case.
Neutral pH Enzymatic Detergent - Low Foam	FHC-404LF-5	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation on instruments and scope washers. 5 gallon plastic jug with handle.
Neutral pH Enzymatic Detergent - Low Foam	FHC-404LF-15	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation on instruments and scope washers. 15 gallon plastic drum.
Neutral pH Enzymatic Detergent - Low Foam	FHC-404LF-30	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation on instruments and scope washers. 30 gallon plastic drum.
Neutral pH Enzymatic Detergent - Low Foam	FHC-404LF-55	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation on instruments and scope washers. 55 gallon plastic drum.
Neutral pH Enzymatic Detergent - Super Concentrate	FHC-404SCLF-2.5	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 2 x 2.5 gallon bottles per case
Neutral pH Enzymatic Detergent - Super Concentrate	FHC-404SCLF-5	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 5 gallon plastic jug w/ handle.
Neutral pH Enzymatic Detergent - Super Concentrate	FHC-404SCLF-15	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 15 gallon plastic drum.
Neutral pH Enzymatic Detergent - Super Concentrate	FHC-404SCLF-30	Surg-ENZ® reduces cleaning time by rapid, efficient removal of stubborn organic encrustation for manual instruments and scope cleaning. 30 gallon plastic drum.

Detergents & Cleaners



Neutral pH Enzymatic Detergent (continued)

Name	Part #	Description
Neutral pH Enzymatic Detergent - Foaming	FHC-405C-1	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 4 x 1 gallon bottle packed 4 per case.
Neutral pH Enzymatic Detergent - Foaming	FHC-405C-5	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 5 gallon plastic jug w/ handle.
Neutral pH Enzymatic Detergent - Foaming	FHC-405C-15	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 15 gallon plastic drum.
Neutral pH Enzymatic Detergent - Foaming	FHC-405C-30	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 30 gallon plastic drum.
Neutral pH Enzymatic Detergent - Foaming	FHC-405C-55	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 55 gallon plastic drum.

Detergents & Cleaners



Neutral pH Enzymatic Detergent (continued)

Name	Part #	Description
Neutral pH Enzymatic Detergent - Low Foam	FHC-405LF-1	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 4 x 1 gallon bottle packed 4 per case.
Neutral pH Enzymatic Detergent - Low Foam	FHC-405LF-5	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 5 gallon plastic jug w/ handle.
Neutral pH Enzymatic Detergent - Low Foam	FHC-405LF-15	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 15 gallon plastic drum.
Neutral pH Enzymatic Detergent - Low Foam	FHC-405LF-30	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 30 gallon plastic drum.
Neutral pH Enzymatic Detergent - Low Foam	FHC-405LF-55	Scope-ZYME® multi-enzymatic cleaner is highly effective for fast, efficient removal of dried blood, proteins, fats, etc.. It is safe for all flexible/rigid scopes, metal, etc. 55 gallon plastic drum.
Enzymatic Scope Cleaning Sponge Set	FHC-442	Scope cleaning sponge set with ready to use sponge soaked in enzymatic detergent for surgical instruments and all scopes.

Detergents & Cleaners



Neutral pH Cart Wash Detergent

Name	Part #	Description
Neutral pH Cart Wash Detergent - Low Foam	FHC-500-5	Cart wash detergent cleans wide range of metal surfaces and is free of metal attacking agents. 5 gallon plastic jug with handle.
Neutral pH Cart Wash Detergent - Low Foam	FHC-500-15	Cart wash detergent cleans wide range of metal surfaces and is free of metal attacking agents. 15 gallon plastic drum.
Neutral pH Cart Wash Detergent - Low Foam	FHC-500-30	Cart wash detergent cleans wide range of metal surfaces and is free of metal attacking agents. 55 gallon plastic drum.
Neutral pH Cart Wash Detergent - Low Foam	FHC-500-55	Cart wash detergent cleans wide range of metal surfaces and is free of metal attacking agents. 55 gallon plastic drum.

Cart Rinse and Drying Concentrate

Name	Part #	Description
Cart Rinse & Drying Agent	FHC-600-5	Cart rinse and drying agent for all brands of cart washers. Designed to provide efficient drying of any cart in any cart washer. 5 gallon plastic jug w/ handle.
Cart Rinse & Drying Agent	FHC-600-15	Cart rinse and drying agent for all brands of cart washers. Designed to provide efficient drying of any cart in any cart washer. 15 gallon plastic drum.
Cart Rinse & Drying Agent	FHC-600-30	Cart rinse and drying agent for all brands of cart washers. Designed to provide efficient drying of any cart in any cart washer. 30 gallon plastic drum.
Cart Rinse & Drying Agent	FHC600-55	Cart rinse and drying agent for all brands of cart washers. Designed to provide efficient drying of any cart in any cart washer. 30 gallon plastic drum.

Detergents & Cleaners



Acidic Stain & Rust Remover

Name	Part #	Description
Stain & Rust Remover	FHC-800-1	Medi-Sheen™ stainless steel reconditioner for manual and machine applications. 4 x 1 gallon bottle packed 4 per case.
Stain & Rust Remover	FHC-800-5	Medi-Sheen™ stainless steel reconditioner for manual and machine applications. 5 gallon plastic jug.
Stain & Rust Remover	FHC-800-15	Medi-Sheen™ stainless steel reconditioner for manual and machine applications. 15 gallon plastic drum.
Stain & Rust Remover	FHC-800-30	Medi-Sheen™ stainless steel reconditioner for manual and machine applications. 30 gallon plastic drum.
Stain & Rust Remover	FHC-800-55	Medi-Sheen™ stainless steel reconditioner for manual and machine applications. 55 gallon plastic drum.
Stain & Rust Remover Spray - Foaming	FHC-800-1S	Medi-Sheen™ stainless steel foaming reconditioner for manual and machine applications. 12 x 22 fl. oz. spray bottles per case (sprayers included)